How to use bimaps from the ".db" annotation packages

Marc Carlson, Herve Pages, Seth Falcon, Nianhua Li

January 29, 2017

1 Introduction

1.0.1 Purpose

AnnotationDbi is used primarily to create mapping objects that allow easy access from R to underlying annotation databases. As such, it acts as the R interface for all the standard annotation packages. Underlying each AnnotationDbi supported annotation package is at least one (and often two) annotation databases. AnnotationDbi also provides schemas for theses databases. For each supported model organism, a standard gene centric database is maintained from public sources and is packaged up as an appropriate organism or "org" package.

1.0.2 Database Schemas

For developers, a lot of the benefits of having the information loaded into a real database will require some knowledge about the database schema. For this reason the schemas that were used in the creation of each database type are included in AnnotationDbi. The currently supported schemas are listed in the DBschemas directory of AnnotationDbi. But it is also possible to simply print out the schema that a package is currently using by using its "_dbschema" method.

There is one schema/database in each kind of package. These schemas specify which tables and indices will be present for each package of that type. The schema that a particular package is using is also listed when you type the name of the package as a function to obtain quality control information.

The code to make most kinds of the new database packages is also included in AnnotationDbi. Please see the vignette on SQLForge for more details on how to make additional database packages.

1.0.3 Internal schema Design of org packages

The current design of the organism packages is deliberately simple and gene centric. Each table in the database contains a unique kind of information and also an internal identifier called _id. The internal _id has no meaning outside of the context of a single database. But _id does connect all the data within a single database.

As an example if we wanted to connect the values in the genes table with the values in the kegg table, we could simply join the two tables using the internal _id column. It is very important to note however that _id does not have any absolute significance. That is, it has no meaning outside of the context of the database where it is used. It is tempting to think that an _id could have such significance because within a single database,

it looks and behaves similarly to an entrez gene ID. But _id is definitely NOT an entrez gene ID. The entrez gene IDs are in another table entirely, and can be connected to using the internal _id just like all the other meaningful information inside these databases. Each organism package is centered around one type of gene identifier. This identifier is found as the gene_id field in the genes table and is both the central ID for the database as well as the foreign key that chip packages should join to.

The chip packages are 'lightweight', and only contain information about the basic probe to gene mapping. You might wonder how such packages can provide access to all the other information that they do. This is possible because all the other data provided by chip packages comes from joins that are performed by AnnotationDbi behind the scenes at run time. All chip packages have a dependency on at least one organism package. The name of the organism package being depended on can be found by looking at its "ORGPKG" value. To learn about the schema from the appropriate organism package, you will need to look at the "_dbschema" method for that package. In the case of the chip packages, the gene_id that in these packages is mapped to the probe_ids, is used as a foreign key to the appropriate organism package.

Specialized packages like the packages for GO and KEGG, will have their own schemas but will also adhere to the use of an internal _id for joins between their tables. As with the organism packages, this _id is not suitable for use as a foreign key.

For a complete listing of the different schemas used by various packages, users can use the available.dbschemas function. This list will also tell you which model organisms are supported.

```
library(DBI)
library(org.Hs.eg.db)
## Loading required package: AnnotationDbi
## Loading required package: stats4
## Loading required package: BiocGenerics
## Loading required package: parallel
##
## Attaching package: 'BiocGenerics'
## The following objects are masked from 'package:parallel':
##
      clusterApply, clusterApplyLB, clusterCall, clusterEvalQ, clusterExport,
##
##
      clusterMap, parApply, parCapply, parLapply, parLapplyLB, parRapply,
##
      parSapply, parSapplyLB
## The following objects are masked from 'package:stats':
##
##
      IQR, mad, xtabs
## The following objects are masked from 'package:base':
##
      Filter, Find, Map, Position, Reduce, any Duplicated, append, as.data.frame,
##
      cbind, colnames, do.call, duplicated, eval, evalq, get, grep, grepl,
##
##
      intersect, is.unsorted, lapply, lengths, mapply, match, mget, order, paste,
##
      pmax, pmax.int, pmin, pmin.int, rank, rbind, rownames, sapply, setdiff,
##
      sort, table, tapply, union, unique, unsplit, which, which.max, which.min
```

```
## Loading required package:
                              Biobase
## Welcome to Bioconductor
##
##
      Vignettes contain introductory material; view with 'browseVignettes()'.
      cite Bioconductor, see 'citation("Biobase")', and for packages
##
      'citation("pkqname")'.
##
## Loading required package:
                              IRanges
## Loading required package:
                              S4Vectors
##
## Attaching package: 'S4Vectors'
## The following objects are masked from 'package:base':
##
##
      colMeans, colSums, expand.grid, rowMeans, rowSums
##
library(AnnotationForge)
available.dbschemas()
```

2 Examples

2.0.1 Basic information

The AnnotationDbi package provides an interface to SQLite-based annotation packages. Each SQLite-based annotation package (identified by a ".db" suffix in the package name) contains a number of AnnDbBimap objects in place of the environment objects found in the old-style environment-based annotation packages. The API provided by AnnotationDbi allows you to treat the AnnDbBimap objects like environment instances. For example, the functions [[, get, mget, and ls all behave the same as they did with the older environment based annotation packages. In addition, new methods like [, toTable, subset and others provide some additional flexibility in accessing the annotation data.

```
library(hgu95av2.db)
##
```

The same basic set of objects is provided with the db packages:

```
ls("package:hgu95av2.db")
    [1] "hgu95av2"
                                 "hgu95av2.db"
##
    [3] "hgu95av2ACCNUM"
                                 "hgu95av2ALIAS2PR0BE"
##
##
    [5] "hgu95av2CHR"
                                 "hgu95av2CHRLENGTHS"
    [7] "hgu95av2CHRLOC"
                                 "hgu95av2CHRLOCEND"
##
    [9] "hgu95av2ENSEMBL"
                                 "hgu95av2ENSEMBL2PR0BE"
## [11] "hgu95av2ENTREZID"
                                 "hgu95av2ENZYME"
## [13] "hgu95av2ENZYME2PR0BE"
                                 "hgu95av2GENENAME"
```

```
## [15] "hgu95av2GO"
                                 "hgu95av2G02ALLPR0BES"
## [17] "hgu95av2G02PR0BE"
                                 "hgu95av2MAP"
## [19] "hgu95av2MAPCOUNTS"
                                 "hgu95av20MIM"
## [21] "hgu95av20RGANISM"
                                 "hgu95av20RGPKG"
## [23] "hgu95av2PATH"
                                 "hgu95av2PATH2PR0BE"
## [25] "hgu95av2PFAM"
                                 "hgu95av2PMID"
## [27] "hgu95av2PMID2PROBE"
                                 "hgu95av2PROSITE"
## [29] "hgu95av2REFSEQ"
                                 "hgu95av2SYMBOL"
## [31] "hgu95av2UNIGENE"
                                 "hgu95av2UNIPROT"
                                 "hgu95av2_dbconn"
## [33] "hgu95av2_dbInfo"
## [35] "hgu95av2_dbfile"
                                 "hgu95av2_dbschema"
```

Exercise 1

Start an R session and use the library function to load the hgu95av2.db software package. Use search() to see that an organism package was also loaded and then use the approriate "_dbschema" methods to the schema for the hgu95av2.db and org.Hs.eg.db packages.

It is possible to call the package name as a function to get some QC information about it.

```
qcdata = capture.output(hgu95av2())
head(qcdata, 20)
    [1] "Quality control information for hgu95av2:"
##
    [2] ""
##
   [3] ""
##
   [4] "This package has the following mappings:"
##
##
   [5] ""
    [6] "hgu95av2ACCNUM has 12625 mapped keys (of 12625 keys)"
##
    [7] "hgu95av2ALIAS2PROBE has 34238 mapped keys (of 120040 keys)"
##
##
    [8] "hgu95av2CHR has 11472 mapped keys (of 12625 keys)"
    [9] "hgu95av2CHRLENGTHS has 93 mapped keys (of 455 keys)"
## [10] "hgu95av2CHRLOC has 11423 mapped keys (of 12625 keys)"
## [11] "hgu95av2CHRLOCEND has 11423 mapped keys (of 12625 keys)"
## [12] "hgu95av2ENSEMBL has 11365 mapped keys (of 12625 keys)"
## [13] "hgu95av2ENSEMBL2PROBE has 9545 mapped keys (of 28015 keys)"
## [14] "hgu95av2ENTREZID has 11474 mapped keys (of 12625 keys)"
## [15] "hgu95av2ENZYME has 2097 mapped keys (of 12625 keys)"
## [16] "hgu95av2ENZYME2PROBE has 779 mapped keys (of 975 keys)"
## [17] "hgu95av2GENENAME has 11474 mapped keys (of 12625 keys)"
## [18] "hgu95av2GO has 11229 mapped keys (of 12625 keys)"
## [19] "hgu95av2GO2ALLPROBES has 18521 mapped keys (of 21161 keys)"
## [20] "hgu95av2G02PROBE has 13971 mapped keys (of 16581 keys)"
```

Alternatively, you can get similar information on how many items are in each of the provided maps by looking at the MAPCOUNTs:

```
hgu95av2MAPCOUNTS
```

To demonstrate the environment API, we'll start with a random sample of probe set IDs.

```
all_probes <- ls(hgu95av2ENTREZID)
length(all_probes)

## [1] 12625

set.seed(0xa1beef)
probes <- sample(all_probes, 5)
probes

## [1] "31882_at" "38780_at" "37033_s_at" "1702_at" "31610_at"</pre>
```

The usual ways of accessing annotation data are also available.

```
hgu95av2ENTREZID[[probes[1]]]

## [1] "9136"

hgu95av2ENTREZID$"31882_at"

## [1] "9136"

syms <- unlist(mget(probes, hgu95av2SYMBOL)))

syms

## 31882_at 38780_at 37033_s_at 1702_at 31610_at

## "RRP9" "AKR1A1" "GPX1" "IL2RA" "PDZK1IP1"
```

The annotation packages provide a huge variety of information in each package. Some common types of information include gene symbols (SYMBOL), GO terms (GO), KEGG pathway IDs (KEGG), ENSEMBL IDs (ENSEMBL) and chromosome start and stop locations (CHRLOC and CHRLOCEND). Each mapping will have a manual page that you can read to describe the data in the mapping and where it came from.

?hgu95av2CHRLOC

Exercise 2

For the probes in 'probes' above, use the annotation mappings to find the chromosome start locations.

2.0.2 Manipulating Bimap Objects

Many filtering operations on the annotation *Bimap* objects require conversion of the *AnnDbBimap* into a *list*. In general, converting to lists will not be the most efficient way to filter the annotation data when using a SQLite-based package. Compare the following two examples for how you could get the 1st ten elements of the hgu95av2SYMBOL mapping. In the 1st case we have to get the entire mapping into list form, but in the second case we first subset the mapping object itself and this allows us to only convert the ten elements that we care about.

```
system.time(as.list(hgu95av2SYMBOL)[1:10])
## vs:
system.time(as.list(hgu95av2SYMBOL[1:10]))
```

There are many different kinds of *Bimap* objects in AnnotationDbi, but most of them are of class *AnnDbBimap*. All /RclassBimap objects represent data as a set of left and right keys. The typical usage of these mappings

is to search for right keys that match a set of left keys that have been supplied by the user. But sometimes it is also convenient to go in the opposite direction.

The annotation packages provide many reverse maps as objects in the package name space for backwards compatibility, but the reverse mappings of almost any map is also available using revmap. Since the data are stored as tables, no extra disk space is needed to provide reverse mappings.

So now that you know about the revmap function you might try something like this:

```
as.list(revmap(hgu95av2PATH)["00300"])
## $`00300`
## [1] "36132_at" "35870_at"
```

Note that in the case of the PATH map, we don't need to use revmap(x) because hgu95av2.db already provides the PATH2PROBE map:

```
x <- hgu95av2PATH
## except for the name, this is exactly revmap(x)
revx <- hgu95av2PATH2PROBE
revx2 <- revmap(x, objName="PATH2PROBE")
revx2
## PATH2PROBE map for chip hgu95av2 (object of class "ProbeAnnDbBimap")
identical(revx, revx2)
## [1] TRUE
as.list(revx["00300"])
## $`00300`
## [1] "36132_at" "35870_at"</pre>
```

Note that most maps are reversible with revmap, but some (such as the more complex GO mappings), are not. Why is this? Because to reverse a mapping means that there has to be a "value" that will always become the "key" on the newly reversed map. And GO mappings have several distinct possibilities to choose from (GO ID, Evidence code or Ontology). In non-reversible cases like this, AnnotationDbi will usually provide a pre-defined reverse map. That way, you will always know what you are getting when you call revmap

While we are on the subject of GO and GO mappings, there are a series of special methods for GO mappings that can be called to find out details about these IDs. Term,GOID, Ontology, Definition,Synonym, and Secondary are all useful ways of getting additional information about a particular GO ID. For example:

```
Term("GD:0000018")
## Loading required package: GD.db
##
## GD:0000018
## "regulation of DNA recombination"
```

```
Definition("GO:0000018")
##
## "Any process that modulates the frequency, rate or extent of DNA recombination, a DNA metaboli
```

Exercise 3

Given the following set of RefSeq IDs: c("NG_005114","NG_007432","NG_008063"), Find the Entrez Gene IDs that would correspond to those. Then find the GO terms that are associated with those entrez gene IDs. org.Hs.eg.db packages.

2.0.3 The Contents and Structure of Bimap Objects

Sometimes you may want to display or subset elements from an individual map. A *Bimap* interface is available to access the data in table (*data.frame*) format using [and toTable.

```
head(toTable(hgu95av2G0[probes]))
##
     probe_id
                   go_id Evidence Ontology
## 1 1702_at GO:0000165
                               TAS
                                         BP
## 2 1702_at GD:0002437
                               IEA
                                         BP
## 3 1702_at GO:0002664
                                         BP
                               IMP
## 4 1702_at GO:0006915
                               TAS
                                         BP
## 5 1702_at GO:0006924
                                         BP
                               IEA
## 6 1702_at GD:0006954
                               IBA
                                         BP
```

The toTable function will display all of the information in a *Bimap*. This includes both the left and right values along with any other attributes that might be attached to those values. The left and right keys of the *Bimap* can be extracted using Lkeys and Rkeys. If is is necessary to only display information that is directly associated with the left to right links in a *Bimap*, then the links function can be used. The links returns a data frame with one row for each link in the bimap that it is applied to. It only reports the left and right keys along with any attributes that are attached to the edge between these two values.

Note that the order of the cols returned by toTable does not depend on the direction of the map. We refer to it as an 'undirected method':

```
toTable(x)[1:6, ]
     probe_id path_id
## 1 1000_at
                04010
## 2 1000_at
                04012
## 3 1000_at
                04062
## 4 1000_at
                04114
## 5 1000_at
                04150
## 6 1000_at
                04270
toTable(revx)[1:6, ]
    probe_id path_id
## 1 1000_at
                04010
## 2 1000_at
                04012
## 3 1000_at
                04062
```

```
## 4 1000_at 04114
## 5 1000_at 04150
## 6 1000_at 04270
```

Notice however that the Lkeys are always on the left (1st col), the Rkeys always in the 2nd col

For length() and keys(), the result does depend on the direction, hence we refer to these as 'directed methods':

```
length(x)
## [1] 12625
length(revx)
## [1] 229
allProbeSetIds <- keys(x)
allKEGGIds <- keys(revx)</pre>
```

There are more 'undirected' methods listed below:

Notice how they give the same result for x and revmap(x)

You might be tempted to think that Lkeys and Llength will tell you all that you want to know about the left keys. But things are more complex than this, because not all keys are mapped. Often, you will only want to know about the keys that are mapped (ie. the ones that have a corresponding Rkey). To learn this you want to use the mappedkeys or the undirected variants mappedLkeys and mappedRkeys. Similarly, the count.mappedkeys, count.mappedLkeys and count.mappedRkeys methods are very fast ways to determine how many keys are mapped. Accessing keys like this is usually very fast and so it can be a decent strategy to subset the mapping by 1st using the mapped keys that you want to find.

```
## [1] "1000_at" "1001_at" "1002_f_at" "1003_s_at" "1004_at"
## [6] "1005_at" "1006_at" "1008_f_at" "1009_at"

count.mappedLkeys(x) # nb of mapped Lkeys
## [1] 9
```

If you want to find keys that are not mapped to anything, you might want to use isNA.

```
y = hgu95av2ENTREZID[isNA(hgu95av2ENTREZID)]  # usage like is.na()
Lkeys(y)[1:4]
## [1] "1007_s_at" "1047_s_at" "1089_i_at" "108_g_at"
```

Exercise 4

How many probesets do not have a GO mapping for the hgu95av2.db package? How many have no mapping? Find a probeset that has a GO mapping. Now look at the GO mappings for this probeset in table form.

2.0.4 Some specific examples

Lets use what we have learned to get information about the probes that are are not assigned to a chromosome:

```
x <- hgu95av2CHR
Rkeys(x)
## [1] "19" "12" "8" "14" "3" "2" "17" "16" "9" "X" "6" "1" "7"
## [14] "10" "11" "22" "5" "18" "15" "Y" "20" "21" "4" "13" "MT" "Un"
chroms <- Rkeys(x)[23:24]
chroms
## [1] "4" "13"
Rkeys(x) <- chroms</pre>
toTable(x)
##
       probe_id chromosome
## 1
      1029_s_at
## 2
       1036_at
                        4
## 3
        1058_at
                       13
## 4
         1065_at
                        13
## 5
        1115_at
                        4
## 6
       1189_at
                        13
## 7
         1198_at
                        13
## 8
        1219_at
                         4
## 9
       1220_g_at
                         4
## 10
       1249_at
                         4
## 11
        1285_at
                         4
## 12
                         4
        1303_at
## 13
       1325_at
                         4
      1348_s_at
## 14
                        13
```

```
## 15
       1369_s_at
## 16
        1377_at
                           4
## 17
       1378_g_at
                           4
## 18
       1451_s_at
                          13
## 19
        1503_at
                          13
## 20
       1507_s_at
                          4
## 21
       1527_s_at
                          13
## 22
        1528_at
                          13
## 23
         1529_at
                          13
## 24
       1530_g_at
                          13
## 25
        1531_at
                          13
## 26
                          13
       1532_g_at
## 27
       1538_s_at
                           4
        1542_at
## 28
                           4
## 29
       1545_g_at
                          13
## 30
        1567_at
                          13
## 31
       1570_f_at
                          13
## 32
       1571_f_at
                          13
## 33
        1593_at
                          4
## 34
         1597_at
                          13
## 35
        1598_g_at
                          13
## 36
         159_at
                           4
## 37
         1600_at
                           4
## 38
         1604_at
                           4
       1605_g_at
## 39
                           4
## 40
        1616_at
                          13
## 41
         1624_at
                           4
## 42
        1629_s_at
                           4
## 43
        1670_at
                          13
       1672_f_at
                          13
## 44
## 45
        1679_at
                           4
## 46
         1708_at
                           4
## 47
       1709_g_at
                           4
## 48
         170_at
                          13
## 49
        1720_at
                           4
## 50
                           4
       1721_g_at
## 51
         1731_at
                           4
## 52
         1732_at
                           4
## 53
         1819_at
                          13
## 54
       1828_s_at
                           4
                           4
## 55
        1836_at
## 56
       1883_s_at
                           4
## 57
       1888_s_at
                           4
## 58
        1900_at
                          13
## 59
       1905_s_at
                          13
## 60
        1913_at
                          4
## 61
         1914_at
                          13
```

```
## 62
       1931_at
                          13
## 63
       1934_s_at
                           4
## 64
        1943_at
                           4
## 65
        1954_at
                           4
## 66
         1963_at
                          13
## 67
       1964_g_at
                          13
## 68
        1987_at
                          4
## 69
         1988_at
                          4
## 70
         1989_at
                          13
## 71
       1990_g_at
                          13
## 72
       2044_s_at
                          13
## 73
         2062_at
                          4
## 74
                           4
       2092_s_at
## 75
        214_at
                           4
## 76
         215_g_at
                           4
## 77
          252_at
                          13
## 78
         253_g_at
                          13
## 79
          260_at
                           4
## 80
         281_s_at
                           4
                           4
## 81
        31314_at
## 82
        31320_at
                          13
## 83
        31333_at
                           4
                           4
## 84
        31345_at
## 85
        31349_at
                           4
## 86
        31356_at
                           4
## 87 31382_f_at
                           4
## 88
        31404_at
                          13
         31408_at
                          4
## 89
## 90
        31464_at
                          13
## 91 31465_g_at
                          13
## 92 31516_f_at
                          13
## 93
        31543_at
                          4
## 94
        31562_at
                          13
## 95
        31584_at
                          13
## 96
      31628_at
                          13
                          4
## 97
      31631_f_at
## 98 31639_f_at
                          13
## 99 31640_r_at
                          13
## 100 31670_s_at
                           4
## 101
                           4
       31684_at
                           4
## 102
       31706_at
## 103 31744_at
                           4
       31753_at
## 104
                          13
       31790_at
                          13
## 105
## 106
        31792_at
                           4
## 107
         31805_at
                           4
## 108 31811_r_at
                           4
```

```
## 109 31847_at
                          13
## 110 31849_at
                          13
## 111
      31851_at
                          13
## 112 31876_r_at
                           4
## 113
        31894_at
                           4
## 114 31969_i_at
                           4
## 115 31970_r_at
                           4
## 116 32006_r_at
## 117 32026_s_at
                           4
## 118
        32080_at
                           4
## 119 32102_at
                          13
## 120
       32145_at
                          4
                           4
## 121 32146_s_at
## 122
        32147_at
                          13
         32148_at
## 123
                          13
## 124 32163_f_at
                           4
## 125 32180_s_at
                          4
## 126
        32220_at
                          13
## 127
        32299_at
                           4
                           4
## 128 32349_at
## 129
        32353_at
                           4
## 130
       32357_at
                           4
## 131
        32368_at
                          13
## 132 32393_s_at
                          4
## 133
        32439_at
                          13
## 134
        32446_at
                          4
## 135
                           4
        32449_at
## 136
        32465_at
                           4
## 137
        32482_at
                          13
## 138
        32506_at
                           4
## 139
        32507_at
                           4
## 140
        32570_at
                           4
## 141
        32580_at
                           4
## 142
        32595_at
                           4
## 143
        32602_at
                           4
## 144
        32641_at
                          13
## 145
        32675_at
                           4
## 146
        32703_at
                           4
## 147
        32768_at
                          13
## 148
                           4
        32769_at
                           4
## 149
        32770_at
## 150
        32771_at
                           4
                           4
## 151
        32812_at
        32822_at
                           4
## 152
## 153
        32832_at
                           4
## 154
        32862_at
                          13
## 155
        32906_at
                          13
```

```
## 156
       32979_at
                           4
                          13
## 157 32986_s_at
## 158
        32998_at
                           4
## 159
        33013_at
                           4
## 160 33068_f_at
                           4
## 161 33069_f_at
                           4
## 162
         33100_at
                           4
## 163
         33150_at
                           4
## 164 33151_s_at
                           4
## 165
        33155_at
                           4
## 166
       33156_at
                           4
## 167
       33168_at
                          13
## 168 33171_s_at
                           4
## 169
         33172_at
                           4
## 170 33173_g_at
                           4
## 171
        33199_at
                          13
## 172
       33208_at
                          13
## 173
       33241_at
                           4
## 174
       33249_at
                           4
## 175
        33267_at
                           4
## 176
        33276_at
                          13
## 177
        33299_at
                          4
## 178
        33318_at
                          13
## 179
        33356_at
                           4
## 180
                           4
        33359_at
## 181
         33369_at
                           4
## 182 33370_r_at
                           4
         33382_at
                           4
## 183
## 184
                           4
        33483_at
## 185
        33488_at
                           4
## 186
        33490_at
                           4
## 187
        33494_at
                           4
## 188
        33519_at
                           4
## 189
        33520_at
                          13
## 190
        33525_at
                           4
                           4
## 191
        33526_at
## 192
                           4
        33529_at
## 193
        33536_at
## 194
        33544_at
                           4
                           4
## 195
        33564_at
## 196
        33576_at
                          13
## 197
        33584_at
                           4
                           4
## 198
        33596_at
## 199
        33657_at
                           4
## 200 33672_f_at
                           4
## 201 33673_r_at
                           4
## 202
         33687_at
                          13
```

```
## 203
         33700_at
                           13
## 204
                           4
         33733_at
## 205
         33791_at
                           13
## 206
         33823_at
                           4
## 207
         33827_at
                           13
## 208
         33837_at
                           4
## 209
         33859_at
                           13
## 210
         33975_at
                           4
## 211
         33990_at
                           4
## 212 33991_g_at
                            4
## 213
         33992_at
                            4
## 214
                           4
         33997_at
## 215
                           4
         34021_at
## 216
         34022_at
                           4
## 217
        34026_at
                           13
## 218
        34029_at
                           4
## 219
        34048_at
                           4
## 220
         34051_at
                           13
## 221
        34058_at
                           4
                            4
## 222
         34075_at
## 223
                            4
         34122_at
## 224
        34131_at
                            4
                            4
## 225
         34144_at
## 226
         34145_at
                            4
## 227
         34149_at
                            4
## 228 34170_s_at
                            4
## 229
                           4
         34181_at
## 230
         34198_at
                           4
## 231
         34211_at
                           13
## 232
         34239_at
                           13
## 233 34240_s_at
                           13
## 234
         34247_at
                           4
## 235
         34248_at
                           4
## 236 34275_s_at
                           4
## 237
         34284_at
                           13
## 238
        34307_at
                           13
## 239
                           4
        34319_at
## 240
        34324_at
                           13
## 241
         34334_at
                           13
## 242
                           13
         34335_at
                           4
## 243
         34341_at
## 244 34342_s_at
                           4
## 245
                           4
       34353_at
## 246
        34398_at
                           13
## 247
        34411_at
                           4
## 248
         34423_at
                           4
## 249
         34459_at
                           13
```

```
## 250 34476_r_at
## 251
                            4
         34482_at
## 252
        34512_at
                            4
## 253
       34551_at
                            4
## 254
        34564_at
                           4
## 255
        34565_at
                           4
## 256
         34578_at
                           13
## 257
         34583_at
                           13
## 258
         34596_at
                           4
## 259 34637_f_at
                           4
## 260 34638_r_at
                           4
## 261
         34657_at
                           13
## 262
         34672_at
                           13
## 263
       34745_at
                           4
## 264
       34803_at
                           13
## 265
         34898_at
                           4
## 266 34953_i_at
                           4
## 267 34954_r_at
                           4
## 268
         34955_at
                           13
                            4
## 269
         34973_at
                           4
## 270
         34984_at
## 271
        34988_at
                            4
                            4
## 272
         35020_at
## 273
         35021_at
                            4
## 274
                            4
         35025_at
## 275
         35028_at
                            4
## 276
                            4
         35039_at
## 277
                            4
         35053_at
## 278
         35061_at
                           4
## 279
         35063_at
                           4
## 280
         35081_at
                           13
## 281
         35105_at
                           13
## 282
        35107_at
                           13
## 283
         35110_at
                           13
## 284
                           4
         35131_at
## 285
                           4
         35134_at
## 286
                           13
         35140_at
## 287
         35147_at
                           13
## 288
         35164_at
                           4
## 289
                           4
         35181_at
                           4
## 290 35182_f_at
## 291
         35193_at
                           13
## 292
                           13
         35213_at
## 293
         35214_at
                           4
## 294
         35215_at
                           4
## 295
         35220_at
                            4
## 296
                            4
         35285_at
```

##	297	35306_at	4
##	298	35344_at	13
##	299	35356_at	4
##	300	35357_at	4
##	301	35371_at	4
##	302	35372_r_at	4
##	303	35400_at	13
##	304	35410_at	4
##	305	35435_s_at	4
##	306	35437_at	4
##	307	35469_at	13
##	308	35470_at	13
##	309	35471_g_at	13
##	310	35481_at	13
##	311	35507_at	4
##	312	35523_at	4
##	313	35554_f_at	13
##	314	35555_r_at	13
##	315	35564_at	4
##	316	35591_at	4
##	317	35656_at	13
##	318	35662_at	4
##	319	35664_at	4
##	320	35678_at	4
##	321	35698_at	4
##	322	35725_at	13
##	323	35730_at	4
##	324	35777_at	4
##	325	35793_at	4
##	326	35827_at	4
##	327	35837_at	4
	328	35845_at	4
##	329	35871_s_at	4
	330	35877_at	13
##	331	35904_at	13
		35939_s_at	13
	333	35940_at	13
	334	35949_at	13
##	335	35972_at	13
	336	35989_at	4
	337	35991_at	4
	338	36012_at	13
	339	36013_at	4
	340	36017_at	13
	341	36021_at	4
	342	36031_at	13
	343	36046_at	4

##	344	36047_at	4
##	345	36065_at	4
##	346	36080_at	4
##	347	36143_at	4
	348	36157_at	4
	349	36188_at	13
	350	36194_at	4
	351	36212_at	13
	352	36243_at	4
		36247_f_at	4
##	354	36269_at	4
##	355	36274_at	13
##	356	36358_at	4
##	357	36363_at	4
##	358	36433_at	4
		36434_r_at	4
	360	36510_at	13
	361	36521_at	13
	362	36606_at	4
	363	36622_at	4
	364	36627_at	4
	365	36659_at	13
##	366	36717_at	4
##	367	36788_at	13
##	368	367_at	13
##	369	36814_at	4
##	370	36830_at	13
	371	36913_at	4
	372	36914_at	4
	373	36915_at	4
	374	36918_at	4
	375	36939_at	4
		36968_s_at	13
	377		4
	378	37006_at	4
##	379	37019_at	4
##	380	37023_at	13
##	381	37056_at	4
	382	37058_at	4
	383	37062_at	4
	384	37067_at	13
	385	37079_at	13
	386	37099_at	13
	387	37109_at	13
	388	37154_at	13
	389	37170_at	4
##	390	37172_at	13

```
## 391
        37173_at
## 392
        37187_at
                           4
## 393
        37206_at
                           4
## 394
        37219_at
                           4
## 395
        37223_at
                           4
## 396
        37243_at
                           4
## 397
         37244_at
                          13
## 398
        37280_at
                           4
## 399
        37282_at
                           4
## 400 37291_r_at
                           4
## 401
         37303_at
                          13
## 402 37322_s_at
                           4
## 403 37323_r_at
                           4
## 404 37356_r_at
                           4
## 405
        37366_at
                           4
## 406
        37404_at
                           4
## 407
       37416_at
                           4
## 408
       37472_at
                           4
## 409
       37518_at
                          13
## 410
        37520_at
                           4
## 411 37521_s_at
                           4
## 412 37522_r_at
                           4
## 413
        37571_at
                          13
## 414
        37578_at
                           4
## 415
       37593_at
                          13
## 416
       37619_at
                           4
## 417
                          13
        37658_at
## 418 37707_i_at
                           4
## 419 37708_r_at
                           4
## 420
        37723_at
                           4
## 421
        37747_at
                           4
## 422
        37748_at
                           4
## 423
        37752_at
                           4
## 424
                          13
        37757_at
## 425
        37767_at
                           4
## 426
                           4
        37840_at
## 427
                           4
        37852_at
## 428
        37926_at
                          13
## 429
        37930_at
                          13
## 430
                           4
        37964_at
                           4
## 431
        38008_at
## 432
         38016_at
                           4
                           4
## 433
         38024_at
                           4
## 434 38025_r_at
## 435
         38035_at
                          13
## 436
         38065_at
                           4
## 437
         38102_at
                          13
```

```
## 438
         38120_at
## 439
                           4
         38168_at
## 440
         38254_at
                           4
## 441 38304_r_at
                           13
## 442
         38353_at
                           13
## 443
         38375_at
                           13
## 444
         38438_at
                           4
## 445
         38485_at
                           4
## 446 38488_s_at
                           4
## 447
         38489_at
                           4
## 448
        38587_at
                           4
## 449
                           4
       38606_at
## 450
       38615_at
                           13
## 451
         38643_at
                           4
## 452
       38649_at
                          13
## 453
        38714_at
                           4
## 454
         38715_at
                           4
## 455
         38736_at
                            4
## 456 38751_i_at
                            4
                            4
## 457 38752_r_at
## 458
         38767_at
                            4
## 459
         38768_at
                            4
                            4
## 460
         38778_at
## 461
         38821_at
                            4
## 462
                            4
         38825_at
## 463
         38838_at
                           4
## 464
                           4
         38854_at
## 465
         38891_at
                           4
## 466
         38957_at
                          13
## 467
         38972_at
                          13
## 468
        38988_at
                           4
## 469
         39028_at
                           13
## 470
        39032_at
                          13
## 471
        39037_at
                           4
## 472
                           4
        39056_at
                           4
## 473
        39083_at
                           13
## 474
         39131_at
## 475
         39132_at
                           4
## 476 39208_i_at
                           4
                           4
## 477 39209_r_at
## 478
         39256_at
                          13
## 479
         39257_at
                          13
## 480
         39269_at
                           13
## 481 39295_s_at
                           4
## 482
        39333_at
                           13
## 483
       39337_at
                           4
## 484
                            4
         39355_at
```

```
## 485
         39369_at
## 486
                           4
         39380_at
## 487
         39382_at
                           4
## 488 39469_s_at
                          13
## 489
         39475_at
                           4
## 490
         39481_at
                           4
## 491
         39488_at
                          13
## 492 39489_g_at
                          13
## 493
        39535_at
                           4
## 494
                           4
         39536_at
## 495
       39554_at
                           4
## 496
                           4
        39555_at
## 497
                           4
         39576_at
## 498
         39579_at
                          13
## 499
         39600_at
                           4
## 500
         39634_at
                           4
## 501 39662_s_at
                           4
## 502
         39665_at
                           4
## 503
        39680_at
                           4
                           4
## 504
         39690_at
         39698_at
## 505
                           4
## 506
        39734_at
                           4
## 507
         39746_at
                           4
## 508
        39748_at
                          13
## 509 39758_f_at
                          13
## 510
        39777_at
                          13
## 511
                           4
         39786_at
        39847_at
                           4
## 512
## 513
                           4
        39850_at
## 514
        39851_at
                           4
## 515
       39852_at
                          13
## 516
        39878_at
                          13
## 517
        39897_at
                           4
## 518
       39924_at
                          13
## 519
        39929_at
                           4
## 520
                           4
       39960_at
## 521
                          13
         39979_at
## 522
         40018_at
                          13
## 523 40058_s_at
                           4
## 524 40059_r_at
                           4
                           4
## 525 40060_r_at
## 526
        40067_at
                          13
## 527
       40072_at
                          13
## 528
        40082_at
                           4
## 529
       400_at
                          13
## 530
         40114_at
                           4
                           4
## 531
         40121_at
```

	532	40148_at	4
	533		13
		40181_f_at	13
	535	40199_at	4
##	536	40217_s_at	4
##	537	40218_at	4
##	538	40225_at	4
##	539	40226_at	4
##	540	40272_at	4
##	541	40310_at	4
##	542	40312_at	13
##	543	40323_at	4
##	544	40349_at	4
##	545	40354_at	4
##	546	40392_at	13
##	547	40404_s_at	13
	548	40449_at	4
	549	40454_at	4
	550	40456_at	4
	551	40473_at	13
	552	40492_at	4
	553	40530_at	4
	554	40570_at	13
		40576_f_at	4
	556	40633_at	13
	557	40681_at	13
	558	40697_at	4
	559	40710_at	4
	560	40710_at	4
	561	40717_at	4
	562	40727_at 40746_at	4
		40770_f_at	4
		40772_at	4
	565	40773_at	4
	566	40818_at	4
	567	40828_at	13
	568	40839_at	13
	569	40853_at	4
		40880_r_at	4
	571	40893_at	13
	572	408_at	4
		40908_r_at	13
	574	40943_at	4
	575	40970_at	13
	576	40989_at	4
	577	40990_at	4
##	578	40991_at	4

```
## 579 40992_s_at
                          4
## 580 40993_r_at
## 581 41014_s_at
                          4
## 582 41024_f_at
                          4
## 583 41025_r_at
                          4
## 584 41026_f_at
                          4
## 585
        41069_at
                         13
## 586 41071_at
                         4
## 587
      41104_at
                         4
## 588 41118_at
                         13
## 589 41119_f_at
                         13
## 590 41145_at
                         4
                         4
## 591
      41148_at
## 592 41182_at
                         13
## 593 41191_at
                         4
## 594
      41276_at
                         13
## 595
      41277_at
                         13
## 596 41300_s_at
                         13
## 597
       41301_at
                         13
                         4
## 598
        41308_at
## 599 41309_g_at
                         4
## 600
        41317_at
                         13
## 601 41318_g_at
                         13
## 602
        41319_at
                         13
## 603 41376_i_at
                         4
## 604 41377_f_at
                          4
## 605 41391_at
                          4
## 606 41392_at
                          4
## 607 41402_at
                          4
## 608 41434_at
                          4
## 609 41436_at
                         13
## 610 41456_at
                         4
## 611
      41459_at
                         13
## 612
        41470_at
                         4
## 613 41491_s_at
                         13
## 614 41492_r_at
                         13
                         13
## 615
       41493_at
## 616
       41534_at
                          4
## 617
       41555_at
                          4
## 618 41556_s_at
                         4
## 619
        41585_at
                         4
## 620 41667_s_at
                         13
## 621 41668_r_at
                         13
## 622
                          4
       41697_at
## 623 41801_at
                          4
## 624 41806_at
                         4
## 625 41860_at
                         13
```

```
## 626
            431_at
                             4
                             4
## 627
            504_at
## 628
          507_s_at
                             4
                             4
## 629
            579_at
## 630
            618_at
                             4
## 631
                             4
            630_at
## 632
         631_g_at
                             4
## 633
            655_at
                             4
## 634
         690_s_at
                             4
## 635
         692_s_at
                             4
## 636
                             4
         764_s_at
## 637
            820_at
                             4
## 638
            886_at
                             4
## 639
            931_at
                            13
## 640
         936_s_at
                             4
## 641
         948_s_at
                             4
## 642
            963_at
                            13
## 643
            975_at
                             4
## 644
            990_at
                            13
## 645
          991_g_at
                            13
```

To get this in the classic named-list format:

```
z <- as.list(revmap(x)[chroms])
names(z)
## [1] "4" "13"
z[["Y"]]
## NULL</pre>
```

Many of the common methods for accessing *Bimap* objects return things in list format. This can be convenient. But you have to be careful about this if you want to use unlist(). For example the following will return multiple probes for each chromosome:

```
chrs = c("12","6")
mget(chrs, revmap(hgu95av2CHR[1:30]), ifnotfound=NA)
## $`12`
## [1] "1018_at" "1019_g_at" "101_at" "1021_at"
##
## $`6`
## [1] "1026_s_at" "1027_at"
```

But look what happens here if we try to unlist that:

```
unlist(mget(chrs, revmap(hgu95av2CHR[1:30]), ifnotfound=NA))
## 121 122 123 124 61 62
## "1018_at" "1019_g_at" "101_at" "1021_at" "1026_s_at" "1027_at"
```

Yuck! One trick that will sometimes help is to use Rfunctionunlist2. But be careful here too. Depending on

what step comes next, Rfunctionunlist2 may not really help you...

```
unlist2(mget(chrs, revmap(hgu95av2CHR[1:30]), ifnotfound=NA))
## 12 12 12 12 6 6
## "1018_at" "1019_g_at" "101_at" "1021_at" "1026_s_at" "1027_at"
```

Lets ask if the probes in 'pbids' mapped to cytogenetic location "18q11.2"?

To coerce this map to a named vector:

The coercion of the reverse map works too but issues a warning because of the duplicated names for the reasons stated above:

```
cyto2pb <- as.character(revmap(x))
## Warning in .local(x, ...): returned vector has duplicated names</pre>
```

2.0.5 Accessing probes that map to multiple targets

In many probe packages, some probes are known to map to multiple genes. The reasons for this can be biological as happens in the arabidopsis packages, but usually it is due to the fact that the genome builds that chip platforms were based on were less stable than desired. Thus what may have originally been a probe designed to measure one thing can end up measuring many things. Usually you don't want to use probes like this, because if they manufacturer doesn't know what they map to then their usefullness is definitely suspect. For this reason, by default all chip packages will normally hide such probes in the standard mappings. But sometimes you may want access to the answers that the manufacturer says such a probe will map to. In such cases, you will want to use the toggleProbes method. To use this method, just call it on a standard mapping and copy the result into a new mapping (you cannot alter the original mapping). Then treat the new mapping as you would any other mapping.

```
## How many probes?
dim(hgu95av2ENTREZID)

## [1] 11473 2

## Make a mapping with multiple probes exposed
```

```
multi <- toggleProbes(hgu95av2ENTREZID, "all")
## How many probes?
dim(multi)
## [1] 13433 2</pre>
```

If you then decide that you want to make a mapping that has only multiple mappings or you wish to revert one of your maps back to the default state of only showing the single mappings then you can use toggleProbes to switch back and forth.

```
## Make a mapping with ONLY multiple probes exposed
multiOnly <- toggleProbes(multi, "multiple")
## How many probes?
dim(multiOnly)

## [1] 1960    2

## Then make a mapping with ONLY single mapping probes
singleOnly <- toggleProbes(multiOnly, "single")
## How many probes?
dim(singleOnly)

## [1] 11473    2</pre>
```

Finally, there are also a pair of test methods hasMultiProbes and hasSingleProbes that can be used to see what methods a mapping presently has exposed.

```
## Test the multiOnly mapping
hasMultiProbes(multiOnly)

## [1] TRUE
hasSingleProbes(multiOnly)

## [1] FALSE

## Test the singleOnly mapping
hasMultiProbes(singleOnly)

## [1] FALSE
hasSingleProbes(singleOnly)
```

2.0.6 Using SQL to access things directly

While the mapping objects provide a lot of convenience, sometimes there are definite benefits to writing a simple SQL query. But in order to do this, it is necessary to know a few things. The 1st thing you will need to know is some SQL. Fortunately, it is quite easy to learn enough basic SQL to get stuff out of a database. Here are 4 basic SQL things that you may find handy:

First, you need to know about SELECT statements. A simple example would look something like this:

```
SELECT * FROM genes;
```

Which would select everything from the genes table.

SELECT gene_id FROM genes;

Will select only the gene_id field from the genes table.

Second you need to know about WHERE clauses:

SELECT gene_id,_id FROM genes WHERE gene_id=1;

Will only get records from the genes table where the gene_id is = 1.

Thirdly, you will want to know about an inner join:

SELECT * FROM genes, chromosomes WHERE genes._id=chromosomes._id;

This is only slightly more complicated to understand. Here we want to get all the records that are in both the 'genes' and 'chromosomes' tables, but we only want ones where the '_id' field is identical. This is known as an inner join because we only want the elements that are in both of these tables with respect to '_id'. There are other kinds of joins that are worth learning about, but most of the time, this is all you will need to do.

Finally, it is worthwhile to learn about the AS keyword which is useful for making long queries easier to read. For the previous example, we could have written it this way to save space:

```
SELECT * FROM genes AS g,chromosomes AS c WHERE g._id=c._id;
```

In a simple example like this you might not see a lot of savings from using AS, so lets consider what happens when we want to also specify which fields we want:

SELECT g.gene_id,c.chromosome FROM genes AS g,chromosomes AS c WHERE g._id=c._id;

Now you are most of the way there to being able to query the databases directly. The only other thing you need to know is a little bit about how to access these databases from R. With each package, you will also get a method that will print the schema for its database, you can view this to see what sorts of tables are present etc.

```
org.Hs.eg_dbschema()
```

To access the data in a database, you will need to connect to it. Fortunately, each package will automatically give you a connection object to that database when it loads.

```
org.Hs.eg_dbconn()
```

You can use this connection object like this:

```
query <- "SELECT gene_id FROM genes LIMIT 10;"
result = dbGetQuery(org.Hs.eg_dbconn(), query)
result</pre>
```

Exercise 5

Retrieve the entrez gene ID and chromosome by using a database query. Show how you could do the same thing by using toTable

2.0.7 Combining data from multiple annotation packages at the SQL level

For a more complex example, consider the task of obtaining all gene symbols which are probed on a chip that have at least one GO BP ID annotation with evidence code IMP, IGI, IPI, or IDA. Here is one way to extract this using the environment-based packages:

```
## Obtain SYMBOLS with at least one GO BP
## annotation with evidence IMP, IGI, IPI, or IDA.
system.time({
bpids <- eapply(hgu95av2GO, function(x) {</pre>
    if (length(x) == 1 && is.na(x))
      NA
    else {
        sapply(x, function(z) {
             if (z$Ontology == "BP")
               z$GOID
             else
               NA
             })
})
bpids <- unlist(bpids)</pre>
bpids <- unique(bpids[!is.na(bpids)])</pre>
g2p <- mget(bpids, hgu95av2G02PR0BE)</pre>
wantedp <- lapply(g2p, function(x) {</pre>
    x[names(x) %in% c("IMP", "IGI", "IPI", "IDA")]
})
wantedp <- wantedp[sapply(wantedp, length) > 0]
wantedp <- unique(unlist(wantedp))</pre>
ans <- unlist(mget(wantedp, hgu95av2SYMBOL))</pre>
})
length(ans)
ans[1:10]
```

All of the above code could have been reduced to a single SQL query with the SQLite-based packages. But to put together this query, you would need to look 1st at the schema to know what tables are present:

```
hgu95av2_dbschema()
```

This function will give you an output of all the create table statements that were used to generate the hgu95av2 database. In this case, this is a chip package, so you will also need to see the schema for the organism package that it depends on. To learn what package it depends on, look at the ORGPKG value:

```
hgu95av2ORGPKG
```

Then you can see that schema by looking at its schema method:

```
org.Hs.eg_dbschema()
```

So now we can see that we want to connect the data in the go_bp, and symbol tables from the org.Hs.eg.sqlite database along with the probes data in the hgu95av2.sqlite database. How can we do that?

It turns out that one of the great conveniences of SQLite is that it allows other databases to be 'ATTACHed'. Thus, we can keep our data in many differnt databases, and then 'ATTACH' them to each other in a modular fashion. The databases for a given build have been built together and frozen into a single version specifically to allow this sort of behavoir. To use this feature, the SQLite ATTACH command requires the filename for the database file on your filesystem. Fortunately, R provides a nice system independent way of getting that information. Note that the name of the database is always the same as the name of the package, with the suffix '.sqlite'.:

```
orgDBLoc = system.file("extdata", "org.Hs.eg.sqlite", package="org.Hs.eg.db")
attachSQL = paste("ATTACH '", orgDBLoc, "' AS orgDB;", sep = "")
dbGetQuery(hgu95av2_dbconn(), attachSQL)
## data frame with 0 columns and 0 rows
```

Finally, you can assemble a cross-db sql query and use the helper function as follows. Note that when we want to refer to tables in the attached database, we have to use the 'orgDB' prefix that we specified in the 'ATTACH' query above.:

```
system.time({
SQL <- "SELECT DISTINCT probe_id,symbol FROM probes, orgDB.gene_info AS gi, orgDB.genes
zz <- dbGetQuery(hgu95av2_dbconn(), SQL)
})

## user system elapsed
## 0.176 0.016 0.381

#its a good idea to always DETACH your database when you are finished...
dbGetQuery(hgu95av2_dbconn(), "DETACH orgDB" )

## data frame with 0 columns and 0 rows</pre>
```

Exercise 6

Retrieve the entrez gene ID, chromosome location information and cytoband information by using a single database query.

Exercise 7

Expand on the example in the text above to combine data from the hgu95av2.db and org.Hs.eg.db with the GO.db package so as to include the GO ID, and term definition in the output.

The version number of R and packages loaded for generating the vignette were:

```
## R version 3.3.2 (2016-10-31)
## Platform: x86_64-pc-linux-gnu (64-bit)
## Running under: Ubuntu 16.04.1 LTS
##
## locale:
## [1] LC_CTYPE=en_US.UTF-8
                                   LC_NUMERIC=C
    [3] LC_TIME=en_US.UTF-8
                                   LC_COLLATE=C
    [5] LC_MONETARY=en_US.UTF-8
                                   LC_MESSAGES=en_US.UTF-8
##
    [7] LC_PAPER=en_US.UTF-8
                                   LC_NAME=C
   [9] LC_ADDRESS=C
##
                                   LC_TELEPHONE=C
## [11] LC_MEASUREMENT=en_US.UTF-8 LC_IDENTIFICATION=C
```

```
## attached base packages:
## [1] parallel stats4 stats
                                    graphics grDevices utils
## [7] datasets methods base
##
## other attached packages:
## [1] GO.db_3.4.0
                              hgu95av2.db_3.2.3
## [3] AnnotationForge_1.16.0 org.Hs.eg.db_3.4.0
## [5] AnnotationDbi_1.36.2 IRanges_2.8.1
## [7] S4Vectors_0.12.1 Biobase_2.34.0 ## [9] BiocGenerics_0.20.0 DBI_0.5-1
## [11] knitr_1.15.1
##
## loaded via a namespace (and not attached):
## [1] Rcpp_0.12.9 XML_3.98-1.5
                                      digest_0.6.12 bitops_1.0-6
## [5] magrittr_1.5 evaluate_0.10 RSQLite_1.1-2 highr_0.6
## [9] stringi_1.1.2 BiocStyle_2.2.1 tools_3.3.2 stringr_1.1.0
## [13] RCurl_1.95-4.8 memoise_1.0.0
```